



IMMEDIATE STABILIZING IMPLANT

DUAL STABILIZATION® LINE


ISI One-Piece Implants

- ➔ Immediate-Load/Function Capability, Same Day Temporization and Function in most cases
- ➔ Perio Lock
- ➔ Cortic-O Thread™
- ➔ Enhanced Surface Treatment
- ➔ Aggressive Stabilizing Thread Pattern
- ➔ Passivated and Sold Sterile
- ➔ Placement in One Procedure
- ➔ Auger-like/Bull Nose Tip




ISI DUAL STABILIZATION® LINE OPTIONS


[3.25mm]



[4.0mm]



[5.0mm]



IOT/O-Ball and CB available in a wide range of lengths:

8mm, 10mm, 12mm, 14mm, 16mm



IMMEDIATE STABILIZING IMPLANT

DUAL STABILIZATION® LINE

INSTRUMENTATION



Catalog# #8 HS • Marking Bur #8 - High Speed



Catalog# 0180 • 1.8 mm Pilot Drill



Catalog# DSO180-10,12,14 • Drill Stops for Pilot Drill



Catalog# ISI-P 3.25, 4.0, 5.0 • Paralleling Pin



Catalog# 0062 • 3.25mm Tissue Punch



Catalog# 0063 • 4.0mm Tissue Punch



Catalog# 0065 • 5.0mm Tissue Punch



Catalog# 3.25 CS • 3.25mm Counter Sink



Catalog# 4.0 CS • 4.0mm Counter Sink



Catalog# 5.0 CS • 5.0mm Counter Sink



Catalog# 0280T • 2.8mm Final Drill for 3.25mm - upper



Catalog# 0300T • 3.0mm Final Drill for 3.25mm - lower



Catalog# 0350T • 3.5mm Final Drill for 4.0mm - upper



Catalog# 0370T • 3.7mm Final Drill for 4.0mm - lower



Catalog# 0450T • 4.5mm Final Drill for 5.0mm - upper



Catalog# 0470T • 4.7mm Final Drill for 5.0mm - lower




Catalog# 0010 • Drill Extender



Catalog# RW • Ratchet Wrench



Catalog# 6101 • Torque Lock Wrench



Catalog# TK Thumb Knob Wrench



Catalog# 0091 4.0 & 5.0mm IOT/O-Ball Driver/Extension



Catalog# 0097 Hand-piece Driver/Adapter



Catalog# 0083N and 0083NS (short) 3.0 & 3.25mm C&B Drivers




Catalog# 0075N and 0075NS (short) 4.0 & 5.0mm C&B Drivers



Catalog# 0076 and 0075S (short) 3.0 & 3.25mm IOT/O-Ball Drivers

RESTORATIVE



Margin Analog Tool

Catalog# MAT 3.25, MAT 4.0, MAT 5.0

Margin Analog

Catalog# 3.25 MA, 4.0 MA, 5.0 MA



C&B Analog

Catalog# ISH-MINIA CB, ISIA 4.0 CB, ISIA 5.0 CB



C&B TRIP (Tissue Retraction Impression Pickup)
Catalog# ISHSDC, ISIC 4.0, ISIC 5.0



Acrylic Coping
Catalog# ISH-MINI AC, 4.0 AC, 5.0 AC

ISI O-BALL ATTACHMENT HEAD [IOT KIT]



ISH-MINI-TSI IOT Kit, ISI 4.0 IOT Kit, or ISI 5.0 IOT Kit includes:

- Encapsulator (003 Cap)
- IOT Analog
- Red Processing O-Ring less retentive (003 Red)
- Black Final O-Ring (003 Blk)

IOT OVERDENTURE MODEL



4 on the Floor™ Model