

Convenient, Simple and Versatile.



OCO Biomedical has one, easy-to-use surgical kit for both the ISI Complete One-piece and TSI Two-piece Dual Stabilization™ Implants. The kit is **organized** and **color-coded** to flow with the implant placement procedures in an **easy to follow, step-by-step format** with all components labeled for easy return to the tray after surgery. The kit's size and low profile allow it to fit easily into the autoclave for **sterilization**. There is also an extra compartment underneath the caddy that allows room for additional components and instruments when necessary.

SURGICAL KIT

Dual Stabilization™ Implant Instrumentation System

- Simple step-by-step drill sequences
- One set of drills for placement of both the ISI & TSI implants
- Color-coded kit and instruments correlate to implant diameter

OCO Biomedical Implant Instrumentation System Components:

- #8 FG Surgical Marking Burr, 1.8 mm Pilot drill and Drill Extension: For marking the site of the osteotomy and setting the initial depth and alignment.
- Paralleling Pins (3.25 mm, 4.00 mm and 5.00 mm): Alignment pins shaped like the abutment; used to check alignment of the osteotomy site.
- Crown & Bridge Drills: For modifying the crown and bridge implant head.
- Tissue Punch, Counter Sink, Upper Final Drill and Lower Final Drill Sets: Drill sequence instrumentation for placement of the 3.00 mm ISD, 3.25 mm ISI & TSI, 4.00 mm ISI & TSI and 5.00 mm ISI & TSI implants (one-piece and two-piece).
- Thumb Knob and Hand-piece Drivers: Manual wrench and hand-piece adaptors for ISD and ISI one-piece implant drivers, .050 mm Hex drivers, and TSI two-piece drivers.
- .050mm Hex drivers (long & short): For engaging healing screws or abutment fixation screws.
- Ratchet and Ratchet Torque Combo Wrenches: For ISD and ISI one-piece implant drivers, .050 mm Hex drivers, and TSI two-piece drivers.
- C&B and "O" Ball Drivers: ISD (3.00 mm) and ISI (3.25 mm, 4.00 mm, 5.00 mm) one-piece implant drivers (C&B or "O" Ball, long or short) for implant placement.
- Extension or 4.00 mm & 5.00 mm Ball Driver: For use as the driver for placement of the ISI "O" Ball implant in 4.00 mm and 5.00 mm diameters, or as a driver extension.
- TSI Drivers: For use as the driver for placement of the TSI two-piece implant (long & short for 3.25 mm, 4.00 mm and 5.00 mm diameter TSI implants). Driver includes the .050 Hex for insertion of healing screw.

