



OCO BIOMEDICAL

Product Catalog

2025/2026




Table of Contents

Complete Implant Surgical System	1
Engage™ Surgical System	2
Guided Surgical System	3
DD-Sinus Lift™ Kit	4
SDI Surgical System	5
Final Osteotomy Former Drill Stop Kit	5
SDI Implants	6
3.0 Implants	7
ISI Implants	8
TSI Implants	9
Engage™ Implants	10
SDI Prosthetics	11
3.0 Prosthetics	11
ISI Prosthetics	12
TSI & ERI Prosthetics	13 - 15
Engage™ Prosthetics	16 - 20
Instrumentation	21 - 23
Regenerative Products	24
Restorative Products	25
Aseptico Implant Unit & Handpiece Packages	26
Aseptico Handpieces	27
Aseptico Accessory Products	27

Complete Implant Surgical System



Places All OCO Biomedical Implants



Complete Implant Surgical System

The Complete Implant Surgical System is organized with a color-coded system that follows the implant placement procedures in an intuitive format. The system will place and restore the SDI, 3.0, ISI, TSI, ERI, and Engage™ OCO Biomedical Implants.

Product Name

Complete Implant Surgical System

Part #

COMPLETE SURGICAL KIT

ORDER TODAY: 800-228-0477

Engage™ Compact Kit

Places OCO Biomedical Engage™ Implants



Engage™ Compact Kit

The OCO Biomedical Engage™ compact kit is organized with the OCO color-coded system that follows the implant placement in an easy to follow step by step format.

Product Name	Part #
Engage™ Compact Kit	ENGAGE COMPACT KIT



Torque Wrench (30-90 Ncm)*

Product Name	Part #
Torque Wrench (30-90 Ncm)	6101
*Sold Separately from Engage™ Compact Kit	

Guided Surgical System

Precision, Predictability, & Performance



Guided Surgical System

The guided surgical system is to be used with 3.25, 4.0, & 5.0mm implant diameters.

Product Name	Part #
Guided Surgical System	GUIDED SURGERY KIT

Guided Surgical Sleeve*



Product Name	Part #
1.8mm Guided Surgical Sleeve	1.8 GSS
3.25mm Guided Surgical Sleeve	3.25 GSS
4.0mm Guided Surgical Sleeve	4.0 GSS
5.0mm Guided Surgical Sleeve	5.0 GSS

*Surgical Sleeves for surgical guide fabrication. To be used with the OCO Biomedical Guided Surgery System.

ORDER TODAY: 800-228-0477

DD-Sinus Lift™ Kit



DD-Sinus Lift Kit™

The OCO Biomedical DD-Sinus Lift™ System contains all of the instrumentation needed to perform a sinus lift procedure using a crestal approach. The DD-Sinus Lift™ technique is unique in that it allows for permeation of the alveolar blood supply from the osteotomy into the sinus cavity for enhanced healing and formation of new bone.

Product Name

DD-Sinus Lift™ Kit

Part #

SINUS LIFT KIT

SDI Surgical System



Places OCO Biomedical SDI Implants (2.2, 2.5, 2.9, & 3.0)

The SDI Surgical system is an economical kit designed specifically for our 2.2, 2.5, 2.9, & 3.0 SDI Implants.



Product Name

SDI Implant Surgical System

Part #

SDI SURGICAL KIT

Final Osteotomy Former Drill Stop Kit

The drill stop kit gives you the vertical control you desire

The Final Osteotomy Former Drill Kit gives you the vertical control and safety you desire for each case. When a drill stop is used, it will prevent over-drilling and ensure that you have gone to the full depth of the desired osteotomy. These stops can be used when placing all styles of OCO Biomedical Implants.



Product Name

Final Osteotomy Former Drill Stop Kit


Part #

DS KIT

ORDER TODAY: 800-228-0477



Small Diameter Implant



Diameter	Length	Collar	Platform	Height
2.2mm*	10, 12, 14, 16mm	2mm	3.0mm	OB - 3.8mm CB - 4.5mm
2.5mm*	10, 12, 14, 16mm	2mm	3.0mm	OB - 3.8mm CB - 4.5mm
2.9mm*	8, 10, 12, 14, 16mm	2mm	3.7mm	OB - 3.8mm CB - 5.5mm

*2.2/2.5/2.9mm O-Ball is 1.8mm & uses 002 O-Ring/002 Encapsulator

2.2mm SDI Crown & Bridge Implant*

Product Name	Part #
2.2mm x 10mm	SDI 2210 CB
2.2mm x 12mm	SDI 2212 CB
2.2mm x 14mm	SDI 2214 CB
2.2mm x 16mm	SDI 2216 CB

2.2mm SDI O-Ball Implant*

Product Name	Part #
2.2mm x 10mm	SDI 2210 OB
2.2mm x 12mm	SDI 2212 OB
2.2mm x 14mm	SDI 2214 OB
2.2mm x 16mm	SDI 2216 OB

2.5mm SDI Crown & Bridge Implant*

Product Name	Part #
2.5mm x 10mm	SDI 2510 CB
2.5mm x 12mm	SDI 2512 CB
2.5mm x 14mm	SDI 2514 CB
2.5mm x 16mm	SDI 2516 CB

2.5mm SDI O-Ball Implant*

Product Name	Part #
2.5mm x 10mm	SDI 2510 OB
2.5mm x 12mm	SDI 2512 OB
2.5mm x 14mm	SDI 2514 OB
2.5mm x 16mm	SDI 2516 OB

2.9mm SDI Crown & Bridge Implant*

Product Name	Part #
2.9mm x 8mm	SDI 2908 CB
2.9mm x 10mm	SDI 2910 CB
2.9mm x 12mm	SDI 2912 CB
2.9mm x 14mm	SDI 2914 CB
2.9mm x 16mm	SDI 2916 CB

2.9mm SDI O-Ball Implant*

Product Name	Part #
2.9mm x 8mm	SDI 2908 OB
2.9mm x 10mm	SDI 2910 OB
2.9mm x 12mm	SDI 2912 OB
2.9mm x 14mm	SDI 2914 OB
2.9mm x 16mm	SDI 2916 OB

Go to page 12 for prosthetic options.

3.0mm Crown & Bridge Implant*

Product Name	Part #
3.0mm x 8mm	3008 CB
3.0mm x 10mm	3010 CB
3.0mm x 12mm	3012 CB
3.0mm x 14mm	3014 CB
3.0mm x 16mm	3016 CB

3.0mm O-Ball Implant*

Product Name	Part #
3.0mm x 8mm	3008 OB
3.0mm x 10mm	3010 OB
3.0mm x 12mm	3012 OB
3.0mm x 14mm	3014 OB
3.0mm x 16mm	3016 OB

3.0



Immediate Stabilizing Implant



Diameter	Length	Collar	Platform	Height
3.0mm	8, 10, 12, 14, 16mm	2mm	3.7mm	OB - 4.5mm CB - 5.5mm

Go to page 11 for prosthetic options.

ORDER TODAY: 800-228-0477



Diameter	Length	Collar	Platform	Height
3.25mm	8, 10, 12, 14, 16mm	2mm	3.7mm	OB - 4.5mm CB - 5.5mm
4.0mm	8, 10, 12, 14, 16mm	2mm	4.7mm	OB - 5.0mm CB - 5.5mm
5.0mm	8, 10, 12, 14, 16mm	2mm	5.7mm	OB - 5.0mm CB - 5.5mm

3.25mm ISI Crown & Bridge Implant*

Product Name	Part #
3.25mm x 8mm	ISI 32508 OB
3.25mm x 10mm	ISI 32510 CB
3.25mm x 12mm	ISI 32512 CB
3.25mm x 14mm	ISI 32514 CB
3.25mm x 16mm	ISI 32516 CB

3.25mm ISI O-Ball Implant*

Product Name	Part #
3.25mm x 8mm	ISI 32508 OB
3.25mm x 10mm	ISI 32510 OB
3.25mm x 12mm	ISI 32512 OB
3.25mm x 14mm	ISI 32514 OB
3.25mm x 16mm	ISI 32516 OB

4.0mm ISI Crown & Bridge Implant*

Product Name	Part #
4.0mm x 8mm	ISI 4008 CB
4.0mm x 10mm	ISI 4010 CB
4.0mm x 12mm	ISI 4012 CB
4.0mm x 14mm	ISI 4014 CB
4.0mm x 16mm	ISI 4016 CB

4.0mm ISI O-Ball Implant*

Product Name	Part #
4.0mm x 8mm	ISI 4008 OB
4.0mm x 10mm	ISI 4010 OB
4.0mm x 12mm	ISI 4012 OB
4.0mm x 14mm	ISI 4014 OB
4.0mm x 16mm	ISI 4016 OB

5.0mm ISI Crown & Bridge Implant*

Product Name	Part #
5.0mm x 8mm	ISI 5008 CB
5.0mm x 10mm	ISI 5010 CB
5.0mm x 12mm	ISI 5012 CB
5.0mm x 14mm	ISI 5014 CB
5.0mm x 16mm	ISI 5016 CB

5.0mm ISI O-Ball Implant*

Product Name	Part #
5.0mm x 8mm	ISI 5008 OB
5.0mm x 10mm	ISI 5010 OB
5.0mm x 12mm	ISI 5012 OB
5.0mm x 14mm	ISI 5014 OB
5.0mm x 16mm	ISI 5016 OB


***ISI Crown & Bridge implant includes the following:**
 (1) Analog, (1) TRIP (Tissue Retraction Impression Pickup),
 (2) Acrylic Copings

***ISI O-Ball implant includes the following:**
 (1) Analog, (1) Encapsulator, (1 each) Red/Black/White O-Ring

Go to page 12 for prosthetic options.



Two-Stage Implant & Esthetic Region Implant




TSI

Diameter	Length	Collar	Platform
3.25mm	8, 10, 12, 14, 16mm	2mm	3.7mm
4.0mm	8, 10, 12, 14, 16mm	2mm	4.7mm
5.0mm	8, 10, 12, 14, 16mm	2mm	5.7mm

ERI

Diameter	Length	Collar	Platform
3.25mm	8, 10, 12, 14, 16mm	1mm	3.7mm
4.0mm	8, 10, 12, 14, 16mm	1mm	4.7mm
5.0mm	8, 10, 12, 14, 16mm	1mm	5.7mm



All ERI Implants have the Internal Hex Connection.

3.25mm TSI/ERI

Product Name	Part #
3.25mm x 8mm*	32508 TSI/ERI
3.25mm x 10mm*	32510 TSI/ERI
3.25mm x 12mm*	32512 TSI/ERI
3.25mm x 14mm*	32514 TSI/ERI
3.25mm x 16mm*	32516 TSI/ERI

*1mm height healing screw included in ERI implant packaging

4.0mm TSI/ERI

Product Name	Part #
4.0mm x 8mm*	4008 TSI/ERI
4.0mm x 10mm*	4010 TSI/ERI
4.0mm x 12mm*	4012 TSI/ERI
4.0mm x 14mm*	4014 TSI/ERI
4.0mm x 16mm*	4016 TSI/ERI

*1mm height healing screw included in ERI implant packaging

5.0mm TSI/ERI

Product Name	Part #
5.0mm x 8mm*	5008 TSI/ERI
5.0mm x 10mm*	5010 TSI/ERI
5.0mm x 12mm*	5012 TSI/ERI
5.0mm x 14mm*	5014 TSI/ERI
5.0mm x 16mm*	5016 TSI/ERI

*1mm height healing screw included in ERI implant packaging

Go to page 13 for prosthetic options.

3.25mm Engage™

Product Name	Part #
3.25mm x 8mm	32508 ENGAGE
3.25mm x 10mm	32510 ENGAGE
3.25mm x 12mm	32512 ENGAGE
3.25mm x 14mm	32514 ENGAGE
3.25mm x 16mm	32516 ENGAGE

3.6mm Engage™ / Engage™ Plus

Product Name	Part #
3.6mm x 10mm	3610 ENGAGE/PLUS
3.6mm x 12mm	3612 ENGAGE/PLUS
3.6mm x 14mm	3614 ENGAGE/PLUS
3.6mm x 16mm	3616 ENGAGE/PLUS

4.0mm Engage™ / Engage™ Plus

Product Name	Part #
4.0mm x 8mm	4008 ENGAGE/PLUS
4.0mm x 10mm	4010 ENGAGE/PLUS
4.0mm x 12mm	4012 ENGAGE/PLUS
4.0mm x 14mm	4014 ENGAGE/PLUS
4.0mm x 16mm	4016 ENGAGE/PLUS

5.0mm Engage™ / Engage™ Plus

Product Name	Part #
5.0mm x 8mm*	5008 ENGAGE/PLUS
5.0mm x 10mm*	5010 ENGAGE/PLUS
5.0mm x 12mm*	5012 ENGAGE/PLUS
5.0mm x 14mm*	5014 ENGAGE/PLUS
5.0mm x 16mm*	5016 ENGAGE/PLUS


6.0mm Engage™ / Engage™ Plus

Product Name	Part #
6.0mm x 6mm	6006 ENGAGE/PLUS
6.0mm x 8mm	6008 ENGAGE/PLUS
6.0mm x 10mm	6010 ENGAGE/PLUS
6.0mm x 12mm	6012 ENGAGE/PLUS
6.0mm x 14mm	6014 ENGAGE

Engage™ / Engage™ Plus



Bone Level Implant / Platform Switching Profile




Engage™ Measurements

Diameter	Length	Platform
3.25mm	8, 10, 12, 14, 16mm	3.0mm
3.6mm	10, 12, 14, 16mm	3.5mm
4.0mm	8, 10, 12, 14, 16mm	3.5mm
5.0mm	8, 10, 12, 14, 16mm	4.5mm
6.0mm	6, 8, 10, 12, 14mm	5.7mm

Engage™ Plus Measurements

Diameter	Length	Platform
4.0mm	8, 10, 12, 14, 16mm	3.5mm
3.6mm	10, 12, 14, 16mm	3.5mm
5.0mm	8, 10, 12, 14, 16mm	3.5mm
6.0mm	8, 10, 12mm	3.5mm



All Engage™ Implants have the Internal Hex Connection.

Go to page 16 for prosthetic options.

ORDER TODAY: 800-228-0477

SDI Prosthetics



Prosthetic Options for SDI



SDI Crown & Bridge/O-Ball Analogs

Product Name	Part #
2.2mm SDI CB/OB Analog	SDIA 2.2/2.5
2.5mm SDI CB/OB Analog	SDIA 2.2/2.5
2.9mm SDI CB/OB Analog	SDIA 2.9



SDI TRIPS (Tissue Retraction Impression Pickups)

Product Name	Part #
2.2mm TRIP	SDIC 2.2/2.5
2.5mm TRIP	SDIC 2.2/2.5
2.9mm TRIP	ISIC 2.9/3/3.25



SDI O-Rings - 10 Packs

Product Name	Part #
002 Red O-Ring	002 RED
002 Black O-Ring	002 BLK
002 White O-Ring	002 WHT
002 Black O-Ring (90D)	002 BN900



SDI Acrylic Copings

Product Name	Part #
2.2mm Acrylic Coping	2.2/2.5 AC
2.5mm Acrylic Coping	2.2/2.5 AC
2.9mm Acrylic Coping	2.9/3.0/3.25 AC



SDI Encapsulators

Product Name	Part #
Encapsulator (002)	002 CAP
Encapsulator (002) w/O-Ring*	002 CAP W/O-RING

*Black O-Ring only



SDI Retaining Rings - 2mm height

Product Name	Part #
2.2mm Retaining Ring	002 RR
2.5mm Retaining Ring	002 RR
2.9mm Retaining Ring	002 RR

3.0 Prosthetics

Prosthetic Options for 3.0



3.0 Crown & Bridge/O-Ball Analogs

Product Name	Part #
3.0mm CB/OB Analog	ISIA 3.0/3.25



3.0 TRIPS (Tissue Retraction Impression Pickups)

Product Name	Part #
3.0mm TRIP	ISIC 2.9/3/3.25



3.0 Retaining Rings - 2mm height

Product Name	Part #
3.0mm Retaining Ring	003 RR



3.0 Acrylic Copings

Product Name	Part #
3.0mm Acrylic Coping	2.9/3.0/3.25 AC



3.0 Encapsulators

Product Name	Part #
Encapsulator (003)	003 CAP
Encapsulator (003) w/O-Ring*	003 CAP W/O-RING

*Black O-Ring only



3.0 O-Rings - 10 Packs

Product Name	Part #
003 Red O-Ring	003 RED
003 Black O-Ring	003 BLK
003 White O-Ring	003 WHT
003 Black O-Ring (90D)	003 BN900



Prosthetic Options for ISI Implants



ISI Crown & Bridge/O-Ball Analogs

Product Name	Part #
3.25mm ISI CB/OB Analog	ISIA 3.0/3.25
4.0mm ISI CB/OB Analog	ISIA 4.0
5.0mm ISI CB/OB Analog	ISIA 5.0



ISI Acrylic Copings

Product Name	Part #
3.25mm ISI Acrylic Coping	2.9/3.0/3.25 AC
3.25mm ISI Acrylic Coping (7mm height)	3.25 AC (7)
4.0mm ISI Acrylic Coping	4.0 AC
4.0mm ISI Acrylic Coping (7mm height)	4.0 AC (7)



ISI TRIPS (Tissue Retraction Impression Pickups)

Product Name	Part #
3.25mm TRIP	ISIC 2.9/3/3.25
4.0mm TRIP	ISIC 4.0
5.0mm TRIP	ISIC 5.0



ISI Encapsulators

Product Name	Part #
Encapsulator (003)	003 CAP
Encapsulator (003) w/O-Ring*	003 CAP W/O-RING

*Black O-Ring only



ISI Retaining Rings - 2mm height

Product Name	Part #
3.0mm Retaining Ring	003 RR



ISI O-Rings - 10 Packs



Product Name	Part #
003 Red O-Ring	003 RED
003 Black O-Ring	003 BLK
003 White O-Ring	003 WHT
003 Black O-Ring (90D)	003 BN900

TSI & ERI Prosthetics



Prosthetic Options for TSI & ERI Implants



TSI/ERI Fixed Solid C & B Abutment - 3.25mm

Product Name	Part #
3.25mm (4.0mm height)	3.25CB 4 TSI
3.25mm (5.5mm height)	3.25CB 5.5 TSI
3.25mm (7mm height)	3.25CB 7 TSI



TSI/ERI Fixed Solid O-Ball Abutment - 3.25mm

Product Name	Part #
3.25mm (0mm height)	325OB0TSI
3.25mm (2mm height)	325OB2TSI
3.25mm (3mm height)	325OB3TSI
3.25mm (4mm height)	325OB4TSI



TSI/ERI Straight Pedestal Abutment

Product Name	Part #
3.25mm diameter	3100 TSI
4.0mm diameter	4100 TSI
5.0mm diameter	4100 TSI



TSI/ERI Straight Shielded Abutment

Product Name	Part #
3.25mm diameter	3150 TSI
4.0mm diameter (short)	4150 TSI S
4.0mm diameter (long)	4150 TSI L
5.0mm diameter (short)	5150 TSI S
5.0mm diameter (long)	5150 TSI L



TSI/ERI Offset Shielded 15° Abutment

Product Name	Part #
3.25mm diameter (short)	3250 TSI S
3.25mm diameter (long)	3250 TSI L
4.0mm diameter (short)	4250 TSI S
4.0mm diameter (long)	4250 TSI L
5.0mm diameter (short)	5250 TSI S
5.0mm diameter (long)	5250 TSI L



TSI/ERI OCO LOK Abutment - 3.25mm

Product Name	Part #
3.25mm (0.7mm height)	1892
3.25mm (2.0mm height)	1893
3.25mm (3.0mm height)	1894
3.25mm (4.0mm height)	1895
3.25mm (5.0mm height)	1896
3.25mm (6.0mm height)	1897



TSI/ERI Fixed Solid C & B Abutment - 4.0/5.0mm

Product Name	Part #
4.0/5.0mm (4mm height)	4CB 4 TSI
4.0/5.0mm (5.5mm height)	4CB 5.5 TSI
4.0/5.0mm (7mm height)	4CB 7 TSI



TSI/ERI Fixed Solid O-Ball Abutment - 4.0/5.0mm

Product Name	Part #
4.0/5.0mm (0mm height)	4OB0TSI
4.0/5.0mm (2mm height)	4OB2TSI
4.0/5.0mm (3mm height)	4OB3TSI
4.0/5.0mm (4mm height)	4OB4TSI



TSI/ERI Offset Pedestal 15° Abutment

Product Name	Part #
3.25mm diameter	3200 TSI
4.0mm diameter	4200 TSI
5.0mm diameter	4200 TSI



TSI/ERI Natural Contour Abutment

Product Name	Part #
3.25mm diameter	3350 TSI-ERI
4.0mm diameter	4350 TSI-ERI
5.0mm diameter	5350 TSI-ERI



TSI/ERI PEEK Temporary Abutments

Product Name	Part #
3.25mm PEEK Abutment	3.25 PA
4.0mm PEEK Abutment	4.0 PA
5.0mm PEEK Abutment	5.0 PA



TSI/ERI OCO LOK Abutment - 4.0mm/5.0mm*

Product Name	Part #
4.0mm (1.0mm height)	8950
4.0mm (2.0mm height)	8951
4.0mm (3.0mm height)	8952
4.0mm (4.0mm height)	8953
4.0mm (5.0mm height)	8954
4.0mm (6.0mm height)	8955

*This is the 4.0 option that works as both 4.0 and 5.0mm options.

TSI & ERI Prosthetics



Prosthetic Options for TSI & ERI



TSI/ERI Healing Abutments - 3.0mm height

Product Name	Part #
3.25mm Healing Abutment	3063 TSI HC
4.0mm Healing Abutment	4063 TSI HC



TSI/ERI Healing Abutments - 5.0mm height

Product Name	Part #
3.25mm Healing Abutment	3065 TSI HC
4.0mm Healing Abutment	4065 TSI HC



TSI/ERI UCLA (Non-Engaging Delrin - Ti Base)

Product Name	Part #
3.25mm TSI/ERI UCLA (non-engaging)	3.25 UCLA NE (TiB)
4.0mm TSI/ERI UCLA (non-engaging)	4.0 UCLA NE (TiB)
5.0mm TSI/ERI UCLA (non-engaging)	5.0 UCLA NE (TiB)



TSI/ERI Delrin UCLA Abutment (Engaging)

Product Name	Part #
3.25mm TSI/ERI UCLA (engaging)	3.25 UCLA E (DEL)
4.0mm TSI/ERI UCLA (engaging)	4.0 UCLA E (DEL)
5.0mm TSI/ERI UCLA (engaging)	5.0 UCLA E (DEL)



TSI/ERI Scannable Abutment

Product Name	Part #
3.25mm TSI/ERI Scannable Abutment	3.25 TiB
4.0mm TSI/ERI Scannable Abutment	4.0 TiB
5.0mm TSI/ERI Scannable Abutment	5.0 TiB



TSI/ERI Indexing Impression Coping*

Product Name	Part #
3.25mm Indexing Impression Coping	3.25 IC
4.0mm Indexing Impression Coping	4.0 IC
5.0mm Indexing Impression Coping	5.0 IC

*Indexing Impression Coping can also be used as Temporary Abutment.



TSI/ERI Healing Abutments - 2.0mm height

Product Name	Part #
3.25mm Healing Abutment	3062 TSI HC
4.0mm Healing Abutment	4062 TSI HC
5.0mm Healing Abutment	5062 TSI HC



TSI/ERI Healing Abutments - 4.0mm height

Product Name	Part #
3.25mm Healing Abutment	3064 TSI HC
4.0mm Healing Abutment	4064 TSI HC
5.0mm Healing Abutment	5064 TSI HC



TSI/ERI UCLA (Engaging Delrin - Ti Base)

Product Name	Part #
3.25mm TSI/ERI UCLA (engaging)	3.25 UCLA E (TiB)
4.0mm TSI/ERI UCLA (engaging)	4.0 UCLA E (TiB)
5.0mm TSI/ERI UCLA (engaging)	5.0 UCLA E (TiB)



TSI/ERI Delrin UCLA Abutment (Non-Engaging)

Product Name	Part #
3.25mm TSI/ERI UCLA (non-engaging)	3.25 UCLA NE (DEL)
4.0mm TSI/ERI UCLA (non-engaging)	4.0 UCLA NE (DEL)
5.0mm TSI/ERI UCLA (non-engaging)	5.0 UCLA NE (DEL)



TSI/ERI Flush Healing Screws

Product Name	Part #
3.25mm	3.25 FLUSH HS
4.0mm	4.0 FLUSH HS
5.0mm	5.0 FLUSH HS



TSI/ERI Implant Analog

Product Name	Part #
3.25mm Implant Analog	TSIA 3.25
4.0mm Implant Analog	TSIA 4.0
5.0mm Implant Analog	TSIA 5.0

ORDER TODAY: 800-228-0477



Prosthetic Options for TSI & ERI Implants



TSI/ERI TRIPS (Tissue Retraction Impression Pickups)

Product Name	Part #
3.25mm TRIP	ISIC 2.9/3/3.25
4.0mm TRIP	ISIC 4.0
5.0mm TRIP	ISIC 5.0



TSI/ERI Offset TRIPS (Tissue Retraction Impression Pickups)

Product Name	Part #
3.25mm Offset TRIP	TSI-OS 3.25
4.0mm Offset TRIP	TSI-OS 4.0
5.0mm Offset TRIP	TSI-OS 5.0



TSI/ERI Retaining Rings - 2mm height

Product Name	Part #
3.0mm Retaining Ring	003 RR



TSI/ERI O-Rings - 10 Packs



Product Name	Part #
003 Red O-Ring	003 RED
003 Black O-Ring	003 BLK
003 White O-Ring	003 WHT
003 Black O-Ring (90D)	003 BN900

OCO-Lok Transfer Pickup - 4 pack

Product Name	Part #
OCO-Lok Transfer Pickup - 4 pack	LOK-PU

OCO-Lok Inserts - 4 pack

Product Name	Part #
Replacement Male, Lab Cap. (Black)	LOK-090
Replacement Male, Extra Soft, 13 lbs, (Yellow)	LOK-088
Replacement Male, Soft, 2.6lbs, (Pink)	LOK-087
Replacement Male, Standard, 40 lbs, (Clear)	LOK-086
Replacement Male, Strong, 5.5 lbs, (Violet)	LOK-089
Metal Housing (2mm height)	LOK-CAP
Protective Disk	LOK-PD



TSI/ERI Solid C& B/O-Ball Analog (5.5mm height)

Product Name	Part #
3.25mm CB/OB Analog	ISIA 3.0/3.25
4.0mm CB/OB Analog	ISIA 4.0



TSI/ERI Acrylic Copings

Product Name	Part #
3.25mm Acrylic Coping	2.9/3.0/3.25 AC
3.25mm Acrylic Coping (7mm height)	3.25 AC (7)
4.0mm Acrylic Coping	4.0 AC
4.0mm Acrylic Coping (7mm height)	4.0 AC (7)
5.0mm Acrylic Coping	5.0 AC



TSI/ERI O-Ball Attachment Kits

Product Name	Part #
3.25mm O-Ball Attachment Kit	3.25 I-OB PK
4.0mm O-Ball Attachment Kit	4.0 TPK-OB
5.0mm O-Ball Attachment Kit	5.0 TPK-OB



OCO-Lok Female Analog - 4 pack

Product Name	Part #
OCO-Lok Analog	LOK-ANALOG



OCO-Lok Processing Kit - 2 pack

Product Name	Part #
OCO-Lok Processing Kit - 2 pack*	8519-2

*Includes denture cap, protective disk, black processing male & replacement males (yellow, pink, clear, purple)



Engage™ Prosthetics



Prosthetic Options for Engage™ Implants



Engage™ Straight Abutment - 3.25mm

Product Name	Part #
3.25mm (1mm collar)	32530-1SA
3.25mm (2mm collar)	32530-2SA
3.25mm (3mm collar)	32530-3SA



Engage™ Straight Abutment - 4.0mm

Product Name	Part #
4.0mm (1mm collar)	4035-1SA
4.0mm (2mm collar)	4035-2SA
4.0mm (3mm collar)	4035-3SA



Engage™ Straight Abutment - 5.0mm

Product Name	Part #
5.0mm (1mm collar)	5045-1SA
5.0mm (2mm collar)	5045-2SA
5.0mm (3mm collar)	5045-3SA



Engage™ Straight Abutment - 6.0mm

Product Name	Part #
6.0mm (1mm collar)	6057-1SA
6.0mm (2mm collar)	6057-2SA
6.0mm (3mm collar)	6057-3SA



Engage™ O-Ball Abutment - 3.25mm

Product Name	Part #
3.25mm (2mm collar)	32530-2OB
3.25mm (3mm collar)	32530-3OB
3.25mm (4mm collar)	32530-4OB



Engage™ O-Ball Abutment - 4.0mm

Product Name	Part #
4.0mm (2mm collar)	4035-2OB
4.0mm (3mm collar)	4035-3OB
4.0mm (4mm collar)	4035-4OB



Engage™ O-Ball Abutment - 5.0mm

Product Name	Part #
5.0mm (2mm collar)	5045-2OB
5.0mm (3mm collar)	5045-3OB
5.0mm (4mm collar)	5045-4OB



Engage™ O-Ball Abutment - 6.0mm

Product Name	Part #
6.0mm (1mm collar)	6057-1OB
6.0mm (2mm collar)	6057-2OB



Engage™ Crown & Bridge Abutment

Product Name	Part #
3.25mm (2mm collar)*	32530-2CB
4.0mm (2mm collar)*	4035-2CB
5.0mm (2mm collar)*	5045-2CB

*For multiple unit restorations only (not for single tooth restoration)



Engage™ Custom Angled 15° Abutment

Product Name	Part #
3.25mm	32530-CAA
4.0mm	4035-CAA
5.0mm	5045-CAA
6.0mm	6057-CAA



Engage™ PEEK Temporary Straight Abutments

Product Name	Part #
3.25mm	32530 PA
4.0mm	4035 PA
5.0mm	5045 PA
6.0mm	6057 PA



Prosthetic Options for Engage™ Implants



Engage™ OCO-LOK Abutments (3.25mm)

Product Name	Part #
3.25mm (0mm height)	1707
3.25mm (1mm height)	1708
3.25mm (2mm height)	1709
3.25mm (3mm height)	1710
3.25mm (4mm height)	1711
3.25mm (5mm height)	1712
3.25mm (6mm height)	1713



Engage™ OCO-LOK Abutments (4.0mm)

Product Name	Part #
4.0mm (0mm height)	8661
4.0mm (1mm height)	8662
4.0mm (2mm height)	8660
4.0mm (2.5mm height)	8663
4.0mm (3.5mm height)	8664
4.0mm (4.5mm height)	8665
4.0mm (5.5mm height)	8625



Engage™ OCO-LOK Abutments (5.0mm)

Product Name	Part #
5.0mm (0mm height)	8671
5.0mm (1mm height)	8672
5.0mm (2.5mm height)	8673
5.0mm (3.5mm height)	8674
5.0mm (4.5mm height)	8675
5.0mm (5.5mm height)	8627
5.0mm (6.5mm height)	8628



Engage™ Healing Abutments Slim (3.0mm height)

Product Name	Part #
3.25mm Engage™ Healing Abutment (Slim)	32530-3HA (Slim)
4.0mm Engage™ Healing Abutment (Slim)	4035-3HA (Slim)
5.0mm Engage™ Healing Abutment (Slim)	5045-3HA (Slim)
6.0mm Engage™ Healing Abutment (Slim)	6057-5HA (Slim)



Engage™ Healing Abutments Slim (5.0mm height)

Product Name	Part #
3.25mm Engage™ Healing Abutment (Slim)	32530-5HA (Slim)
4.0mm Engage™ Healing Abutment (Slim)	4035-5HA (Slim)
5.0mm Engage™ Healing Abutment (Slim)	5045-5HA (Slim)
6.0mm Engage™ Healing Abutment (Slim)	6057-5HA (Slim)



Engage™ Healing Abutments (3.0mm height)

Product Name	Part #
3.25mm Engage™ Healing Abutment	32530-3HA
4.0mm Engage™ Healing Abutment	4035-3HA
5.0mm Engage™ Healing Abutment	5045-3HA
6.0mm Engage™ Healing Abutment	6057-5HA



Engage™ Healing Abutments (5.0mm height)

Product Name	Part #
3.25mm Engage™ Healing Abutment	32530-5HA
4.0mm Engage™ Healing Abutment*	4035-5HA
5.0mm Engage™ Healing Abutment	5045-5HA
6.0mm Engage™ Healing Abutment	6057-5HA



Engage™ Flush Healing Screws

Product Name	Part #
3.25mm Engage™ Flush Healing Screw	32530-HS
4.0mm Engage™ Flush Healing Screw	4035-HS
5.0mm Engage™ Flush Healing Screw	5045-HS
6.0mm Engage™ Flush Healing Screw	6057-HS

Engage™ Prosthetics



Prosthetic Options for Engage™ Implants



Engage™ 4.0mm Multi-Unit Abutment*

Product Name	Part #
Multi-Unit Abutment (1mm height)	4035-1MUA
Multi-Unit Abutment (2mm height)	4035-2MUA

*3.5mm platform



Engage™ 4.0mm Multi-Unit Angled 17° Abutment*

Product Name	Part #
17° Multi-Unit Abut. (2mm height)	4035-2MUA (Angled 17)
17° Multi-Unit Abut. (3mm height)	4035-3MUA (Angled 17)

*3.5mm platform



Engage™ 4.0mm Multi-Unit Angled 30° Abutment*

Product Name	Part #
30° Multi-Unit Abut. (3mm height)	4035-3MUA (Angled 30)

*3.5mm platform



Engage™ 5.0mm Multi-Unit Abutment*

Product Name	Part #
Multi-Unit Abutment (1mm height)	5045-1MUA
Multi-Unit Abutment (2mm height)	5045-2MUA
Multi-Unit Abutment (3mm height)	5045-3MUA
Multi-Unit Abutment (5mm height)	5045-5MUA

*4.5mm platform



Engage™ 5.0mm Multi-Unit Angled 17° Abutment*

Product Name	Part #
17° Multi-Unit Abut. (2mm height)	5045-2MUA (Angled 17)
17° Multi-Unit Abut. (3mm height)	5045-3MUA (Angled 17)
17° Multi-Unit Abut. (5mm height)	5045-5MUA (Angled 17)

*4.5mm platform



Engage™ 5.0mm Multi-Unit Angled 30° Abutment*

Product Name	Part #
30° Multi-Unit Abut. (3mm height)	5045-3MUA (Angled 30)
30° Multi-Unit Abut. (5mm height)	5045-5MUA (Angled 30)

*4.5mm platform



Engage™ Multi-Unit Abutment Analog

Product Name	Part #
Multi-Unit Abutment Analog	MUA ANALOG



Engage™ Multi-Unit Gold Bar Coping

Product Name	Part #
Multi-Unit Gold Bar Coping	MU BAR COPING



Engage™ Multi-Unit Gold/Plastic Coping

Product Name	Part #
Multi-Unit Gold/Plastic Coping	MU COPING



Engage™ Multi-Unit Inclusive Guide Pin

Product Name	Part #
Multi-Unit Inclusive Guide Pin	MU GUIDE PIN



Engage™ Multi-Unit Driver

Product Name	Part #
Multi-Unit Driver	MU DRIVER




Engage™ Multi-Unit Temporary Healing Cap

Product Name	Part #
Multi-Unit Temporary Healing	MU HEALING CAP



Engage™ Multi-Unit Prosthetic Screw

Product Name	Part #
Multi-Unit Prosthetic Screw	MU SCREW



Engage™ Multi-Unit Titanium Temporary Abutment*

Product Name	Part #
Multi-Unit Temporary Healing	MU TEMP

*Includes prosthetic screw



Engage™ Multi-Unit Scanning Abutment*

Product Name	Part #
Multi-Unit Scanning Abutment	MU SCAN

*Includes prosthetic screw



Engage™ Multi-Unit Impression Coping (Open Tray)

Product Name	Part #
Multi-Unit Impression Coping (O)	MU-IC (OPEN TRAY)



Engage™ Multi-Unit Impression Coping (Closed Tray)

Product Name	Part #
Multi-Unit Impression Coping (C)	MU-IC (CLOSED TRAY)

ORDER TODAY: 800-228-0477

18

Engage™ Prosthetics

Prosthetic Options for Engage™ Implants



Engage™ UCLA (Non-Engaging Delrin - Ti Base)*

Product Name	Part #
3.25mm Engage™ UCLA (non-engaging)	32530 UCLA NE (TiB)
4.0mm Engage™ UCLA (non-engaging)	4035 UCLA NE (TiB)
5.0mm Engage™ UCLA (non-engaging)	5045 UCLA NE (TiB)
6.0mm Engage™ UCLA (non-engaging)	6057 UCLA NE (TiB)
7.0mm Engage™ UCLA (non-engaging)	6057 UCLA NE (TiB)
8.0mm Engage™ UCLA (non-engaging)	6057 UCLA NE (TiB)

*Also available with a Gold Base. Please call (800) 228-0477 for additional information.



Engage™ UCLA (Engaging Delrin - Ti Base)*

Product Name	Part #
3.25mm Engage™ UCLA (engaging)	32530 UCLA E (TiB)
4.0mm Engage™ UCLA (engaging)	4035 UCLA E (TiB)
5.0mm Engage™ UCLA (engaging)	5045 UCLA E (TiB)
6.0mm Engage™ UCLA (engaging)	6057 UCLA E (TiB)
7.0mm Engage™ UCLA (engaging)	6057 UCLA E (TiB)
8.0mm Engage™ UCLA (engaging)	6057 UCLA E (TiB)

*Also available with a Gold Base. Please call (800) 228-0477 for additional information.



Engage™ UCLA (Non-Engaging Delrin - Delrin Base)

Product Name	Part #
3.25mm Engage™ UCLA (non-engaging)	32530 UCLA NE (DEL)
4.0mm Engage™ UCLA (non-engaging)	4035 UCLA NE (DEL)
5.0mm Engage™ UCLA (non-engaging)	5045 UCLA NE (DEL)
6.0mm Engage™ UCLA (non-engaging)	6057 UCLA NE (DEL)



Engage™ UCLA (Engaging Delrin - Delrin Base)

Product Name	Part #
4.0mm Engage™ UCLA (engaging)	4035 UCLA E (DEL)
5.0mm Engage™ UCLA (engaging)	5045 UCLA E (DEL)
6.0mm Engage™ UCLA (engaging)	6057 UCLA E (DEL)



Engage™ Digital Scan Body*

Product Name	Part #
3.25mm Engage™ Scan Body	32530 SB
4.0mm Engage™ Scan Body*	4035 SB
5.0mm Engage™ Scan Body	5045 SB
6.0mm Engage™ Scan Body	6057 SB

*Compatible with the Engage™ Plus 4.0, 5.0, and 6.0mm



ENGAGE™ TIBASE (SHORT)

Product Name	Part #
3.25mm	32350 TiBase (SHORT)
4.0mm	4035 TiBase (SHORT)
5.0mm	5045 TiBase (SHORT)
6.0mm	6057 TiBase (SHORT)



ENGAGE™ TIBASE (LONG)

Product Name	Part #
3.25mm	32350 TiBase (LONG)
4.0mm	4035 TiBase (LONG)
5.0mm	5045 TiBase (LONG)
6.0mm	6057 TiBase (LONG)



Engage™ Scannable Abutment

Product Name	Part #
3.25mm Engage™ Scannable Abutment	32530-TiB
4.0mm Engage™ Scannable Abutment	4035-TiB
5.0mm Engage™ Scannable Abutment	5045-TiB
6.0mm Engage™ Scannable Abutment	6057-TiB



Engage™ Impression Coping*

Product Name	Part #
3.25mm Engage™ Impression Coping	32530 IC/TA
4.0mm Engage™ Impression Coping	4035 IC/TA
5.0mm Engage™ Impression Coping	5045 IC/TA
6.0mm Engage™ Impression Coping	6057 IC/TA



Engage™ Digital Analog

Product Name	Part #
3.25mm Engage™ Analog	32530DA
4.0mm Engage™ Analog	4035DA
5.0mm Engage™ Analog	5045DA
6.0mm Engage™ Analog	6057DA



Engage™ Analog

Product Name	Part #
3.25mm Engage™ Analog	32530A
4.0mm Engage™ Analog	4035A
5.0mm Engage™ Analog	5045A
6.0mm Engage™ Analog	6057A



Engage™ Retaining Rings - 2mm height

Product Name	Part #
3.0mm Retaining Ring	003 RR

Engage™ Prosthetics

Prosthetic Options for Engage™ Implants



MILLING BLANK

Product Name	Part #
3.25mm Milling Blank (3.0mm Platform)	32530 MB
4.0mm Milling Blank (3.5mm Platform)*	4035 MB
5.0mm Milling Blank (4.5mm Platform)*	5045 MB
6.0mm Milling Blank (5.7mm Platform)	6057 MB



OCO-Lok Female Analog - 2 pack

Product Name	Part #
OCO-Lok Analog	LOK-ANALOG



OCO-Lok Impression Coping - 2 pack

Product Name	Part #
OCO-Lok Transfer Pickup - 4 pack	LOK-PU



OCO-Lok Processing Kit - 2 pack*

Product Name	Part #
OCO-Lok Processing Kit - 2 pack*	8519-2

*Includes denture cap, protective disk, black processing male & replacement males (yellow, pink, clear, purple)

OCO-Lok Inserts - 4 pack

Product Name	Part #
Replacement Male, Lab Cap. (Black)	LOK-090
Replacement Male, Extra Soft, 1.3 lbs, (Yellow)	LOK-088
Replacement Male, Soft, 2.6lbs, (Pink)	LOK-087
Replacement Male, Standard, 4.0 lbs, (Clear)	LOK-086
Replacement Male, Strong, 5.5 lbs, (Violet)	LOK-089
<hr/>	
Metal Housing (2mm height)	LOK-CAP
Protective Disk	LOK-PD



Pilot Drills

Product Name	Part #
1.5mm Pilot Drill*	0150
1.8mm Pilot Drill	0180
2.0 Pilot Drill 8mm Length	0200-8
2.0 Pilot Drill 10mm Length	0200-10
2.0 Pilot Drill 12mm Length	0200-12
2.0 Pilot Drill 14mm Length	0200-14
2.0 Pilot Drill 16mm Length	0200-16

*Use the 1.5mm Pilot Drill for the 2.2mm SDI



Drill Extension

Product Name	Part #
Drill Extension	0010



#8 Marking Burs

Product Name	Part #
#8 High Speed Marking Bur	#8 HS
#8 Latch Marking Bur	#8 LT



Adjustable Depth Gauge

Product Name	Part #
Adjustable Depth Gauge	ADJ. DEPTH GAUGE



3.25mm Tissue Punch

Product Name	Part #
3.25mm Tissue Punch	0062



4.0mm Tissue Punch

Product Name	Part #
4.0mm Tissue Punch	0063



5.0mm Tissue Punch

Product Name	Part #
5.0mm Tissue Punch	0065



6.0mm Tissue Punch

Product Name	Part #
6.0mm Tissue Punch	0060



Drill Stops for Pilot Drills

Product Name	Part #
6mm Drill Stop	DS0180-06
8mm Drill Stop	DS0180-08
10mm Drill Stop	DS0180-10
12mm Drill Stop	DS0180-12
14mm Drill Stop	DS0180-14
Drill Stop Kit (all 5 drill stops)	PILOT DS KIT



Paralleling Pins

Product Name	Part #
3.25mm Paralleling Pin	ISI-P 3.25
4.0mm Paralleling Pin	ISI-P 4.0
5.0mm Paralleling Pin	ISI-P 5.0
Engage	ENG-PP



Profiling Bur for 3.0/ISI/TSI/ERI

Product Name	Part #
3.0/3.25mm Profiling Bur	0098
4.0mm Profiling Bur	0099
5.0mm Profiling Bur	0100



2.4mm Osteotomy Former/Final Drill*

Product Name	Part #
2.4mm Osteotomy Former	0240T

*2.4mm diameter used for 3.0(mm) implants



2.8mm Osteotomy Former/Final Drill*

Product Name	Part #
2.8mm Osteotomy Former	0280T

*2.8mm diameter used for 3.25mm ISI/TSI/ERI/Engage™ - Maxillary drill



3.0mm Osteotomy Former/Final Drill*

Product Name	Part #
3.0mm Osteotomy Former	0300T

*3.0mm diameter used for 3.25mm ISI/TSI/ERI/Engage™ - Mandibular drill
*3.0mm diameter used for 3.6mm Engage™ Maxillary Drill



3.3mm Osteotomy Former/Final Drill*

Product Name	Part #
3.3mm Osteotomy Former	0330T

*3.3mm diameter used for 3.6mm Engage™ - Mandibular drill



3.5mm Osteotomy Former/Final Drill*

Product Name	Part #
3.5mm Osteotomy Former	0350T

*3.5mm diameter used for 4.0mm ISI/TSI/ERI/Engage™ - Maxillary drill



3.7mm Osteotomy Former/Final Drill*

Product Name	Part #
3.7mm Osteotomy Former	0370T

*3.7mm diameter used for 4.0mm ISI/TSI/ERI/Engage™ - Mandibular drill



4.5mm Osteotomy Former/Final Drill*

Product Name	Part #
4.5mm Osteotomy Former	0450T

*4.5mm diameter used for 5.0mm ISI/TSI/ERI/Engage™ - Maxillary drill



5.6mm Osteotomy Former/Final Drill*

Product Name	Part #
5.6mm Osteotomy Former	0560T

*5.6mm diameter used for 6.0mm Engage™



Crestal Former Drill

Product Name	Part #
3.25mm Crestal Former Drill	32530CF
4.0mm Crestal Former Drill	4035CF
5.0mm Countersink Drill	5045CF
6.0mm Countersink Drill	6057CF



2.2/2.5mm SDI Driver*

Product Name	Part #
2.2/2.5mm SDI Driver (long)	MD-L
2.2/2.5mm SDI Driver (short)	MD-S

*Places 2.2/2.5mm SDI implant both C & B and O-Ball



4.0/5.0mm ISI O-Ball Driver/Extension*

Product Name	Part #
ISI O-Ball Driver/Extension	0091

*Places 4.0/5.0mm ISI O-Ball implants & can also used as an extension for drivers



4.0/5.0mm Crown & Bridge Driver*

Product Name	Part #
4.0/5.0mm C & B Driver (long)	0075N
4.0/5.0mm C & B Driver (short)	0075NS

*Places 4.0/5.0mm ISI C & B implants as well as 4.0/5.0mm C & B abutments for TSI/ERI



Handpiece Driver/Adapter

Product Name	Part #
Handpiece Driver/Adapter	0097



4.7mm Osteotomy Former/Final Drill*

Product Name	Part #
4.7mm Osteotomy Former	0470T

*4.7mm diameter used for 5.0mm ISI/TSI/ERI/Engage™ - Mandibular drill



Countersink Drills*

Product Name	Part #
3.25mm Countersink Drill	3.25 CS
4.0mm Countersink Drill	4.0 CS
5.0mm Countersink Drill	5.0 CS

*For use with OCO Biomedical tissue level implants (ISI/TSI/ERI)



2.9/3.0/3.25mm O-Ball Driver*

Product Name	Part #
2.9/3.0/3.25mm O-Ball Driver (long)	0076
2.9/3.0/3.25mm O-Ball Driver (short)	0076S

*Places 2.9 SDI, 3.0mm, 3.25mm ISI O-Ball implants as well as 3.25mm O-Ball Abutments for 3.25/4.0mm TSI/ERI



2.9/3.0/3.25mm Crown & Bridge Driver*

Product Name	Part #
2.9/3.0/3.25mm C&B Driver (long)	0083N
2.9/3.0/3.25mm C&B Driver (short)	0083NS

*Places 2.9 SDI, 3.0mm, 3.25mm ISI C & B implants as well as 3.25mm C & B abutments for TSI/ERI.



.050mm Hex Driver*

Product Name	Part #
.050mm Hex Driver (long)	.050 L
.050mm Hex Driver (short)	.050 S

*Places healing abutments, locking screws and cover screws



Internal Hex Driver (TSI/ERI/Engage™)

Product Name	Part #
3.25mm Internal Hex Driver (short)	0017
3.25mm Internal Hex Driver (long)	0018
4/5.0mm Internal Hex Driver (short)	0016
4/5.0mm Internal Hex Driver (long)	0015
6.0mm Internal Hex Driver (short)*	0020
6.0mm Internal Hex Driver (long)*	0019



Osteotomy Former Drill Stops

Product Name	Part #
8mm Length (for 2.4, 2.8, 3.0mm*)	DS 2.4/2.8/3.0-08
10mm Length	DS 2.4/2.8/3.0-10
12mm Length	DS 2.4/2.8/3.0-12
14mm Length	DS 2.4/2.8/3.0-14
16mm Length	DS 2.4/2.8/3.0-16
8mm Length (for 3.5, 3.7mm*)	DS 3.5/3.7-08
10mm Length	DS 3.5/3.7-10
12mm Length	DS 3.5/3.7-12
14mm Length	DS 3.5/3.7-14
16mm Length	DS 3.5/3.7-16
8mm Length (for 4.5, 4.7mm*)	DS 4.5/4.7-08
10mm Length	DS 4.5/4.7-10
12mm Length	DS 4.5/4.7-12
14mm Length	DS 4.5/4.7-14
16mm Length	DS 4.5/4.7-16
6mm Length (for 5.6mm*)	DS 5.6-06
8mm Length	DS 5.6-08
10mm Length	DS 5.6-10

*Osteotomy Former diameters



Ratchet Wrench

Product Name	Part #
Ratchet Wrench	RW



Thumb Knob Wrench

Product Name	Part #
Thumb Knob Wrench	TK



Surgical Driver

Product Name	Part #
Surgical Driver	SD



OCO-Lok Insertion Tool

Product Name	Part #
OCO-Lok Insertion Tool	LOK-TOOL



Ridge Mapping Caliper

Product Name	Part #
Ridge Mapping Caliper	RMC



Torque Wrench (30-90 Ncm)

Product Name	Part #
Torque Wrench (30-90 Ncm)	6101



Adjustable Torque Wrench (30-90 Ncm)

Product Name	Part #
Adjustable Torque Wrench	ATW



Thumb Knob Wrench (for Torque Wrench)

Product Name	Part #
Thumb Knob Wrench (for Torque Wrench)	TK (6101)



Abutment Prepping Tool

Product Name	Part #
Abutment Prepping Tool	ABUTMENT TOOL



Offset Osteotomes

Product Name	Part #
1.8mm Offset Osteotome	OST-OS 1.8
2.65mm Offset Osteotome	OST-OS 2.65
3.4mm Offset Osteotome	OST-OS 3.4
4.4mm Offset Osteotome	OST-OS 4.4
OST-OS Kit (all 4 w/case)	OST-OS Kit



Concave Tip Offset Osteotomes

Product Name	Part #
2.65 Concave Tip Offset Osteotome	OST-OS-C 2.65
3.4 Concave Tip Offset Osteotome	OST-OS-C 3.40
4.4 Concave Tip Offset Osteotome	OST-OS-C 4.40
OST-OS-C Kit (all 3 w/case)	OST-OS-C Kit

Regenerative Products

Allograft



Mineralized Cortical/ Cancellous Blend (250-1000)

Product Name	Part #
0.5cc Mineralized Cortico-Cancellous Blend	MCC05
1.0cc Mineralized Cortico-Cancellous Blend	MCC10
2.0cc Mineralized Cortico-Cancellous Blend	MCC20
2.5cc Mineralized Cortico-Cancellous Blend	MCC25
5.0cc Mineralized Cortico-Cancellous Blend	MCC50



Allograft DBM Putty

Product Name	Part #
Allograft DBM Putty .7cc Syringe	DBM07
Allograft DBM Putty 1.0cc Syringe	DBM10
Allograft DBM Putty 2.5cc Syringe	DBM25
Allograft DBM Putty 5.0cc Syringe	DBM50



Allograft DBM Crunch Putty

Product Name	Part #
Allograft DBM Crunch Putty .7cc Syringe	DBM07C
Allograft DBM Crunch Putty 1.0cc Syringe	DBM10C
Allograft DBM Crunch Putty 2.5cc Syringe	DBM25C
Allograft DBM Crunch Putty 5.0cc Syringe	DBM50C

Alloplast

OSTEOGEN® Grafting Plugs

Product Name	Part #
Bone Grafting Plug 6mm X 25 cm (5-Pack)	OP 625
Bone Grafting Plug 10mm X 20 cm (5-Pack)	OP 1020
Bone Grafting Plug 6mm X 25 cm (10-Pack)	OP 625-10
Bone Grafting Plug 10mm X 20 cm (10-pack)	OP 1020-10



Collagen Products



Curagen Plugs: Collagen Plug

Product Name	Part #
Curagen Plugs: Collagen Plug 10X20 (10PK)	HAD001020



InterCollagen® Membrane

Product Name	Part #
Resorbable Collagen Membrane 15 X 20	ICG1520
Resorbable Collagen Membrane 20 X 30	ICG2030
Resorbable Collagen Membrane 30 X 40	ICG3040

*Shipping Costs Vary For Regenerative Products

ORDER TODAY: 800-228-0477

24

WWW.OCOBIOMEDICAL.COM



AIT Kit (002)*

Product Name	Part #
AIT Kit (002)*	002 O-RING KIT

*Kit contains: (1) Straight Post (002 POST), (1) Encapsulator (002 CAP), (2) Red O-Rings (002 RED), (1) Black O-Ring (002 BLK)



AIT Kit (003)*

Product Name	Part #
AIT Kit (003)*	003 O-RING KIT

*Kit contains: (1) Straight Post (003 POST), (1) Encapsulator (003 CAP), (2) Red O-Rings (003 RED), (1) Black O-Ring (003 BLK)



Encapsulators

Product Name	Part #
Encapsulator (002)	002 CAP
Encapsulator (002) w/O-Ring*	002 CAP W/O-RING
Encapsulator (003)	003 CAP
Encapsulator (003) w/O-Ring*	003 CAP W/O-RING
Encapsulator (005)	005 CAP

*Black O-Ring only



Retaining Rings - 2mm height

Product Name	Part #
2.0mm Retaining Ring	002 RR
3.0mm Retaining Ring	003 RR



Straight Posts

Product Name	Part #
002 Straight Post	002 POST
003 Straight Post	003 POST
005 Straight Post	005 POST



O-Rings - 10 Packs

Product Name	Part #
002 Red O-Ring	002 RED
002 Black O-Ring	002 BLK
002 White O-Ring	002 WHT
002 Black O-Ring (90D)	002 BN900
003 Red O-Ring	003 RED
003 Black O-Ring	003 BLK
003 White O-Ring	003 WHT
003 Black O-Ring (90D)	003 BN900



Distal Extension Posts

Product Name	Part #
002 Distal Extension Post	002 DX POST
003 Distal Extension Post	003 DX POST

Implant Unit

Aseptico Product Packages



Aseptico Implant Motor & Surgical Motor Packages (AEU-1000)

Product Name

AEU-1000-OCO



Aseptico Implant Motor & Surgical Motor Packages (AEU-6000)

The AEU-6000 is quickly becoming the best-selling implant motor on the market due to its reliability and affordability.

This versatile implant and surgical motor performs standard or mini implants, third molar extractions and other surgical applications, and even endodontics with auto-stop-reverse.

Product Name

AEU-6000/Mont Blanc Package (Non-fiber optic)

AEU-6000/Impulsion Package



Aseptico Implant Motor & Surgical Motor Packages (AEU-7000)

With more power, features, and versatility than competing systems, the AEU-7000 is the most reliable, accurate, and comprehensive implant, surgical & endo motor on the market today.

Product Name

AEU-7000/Mont Blanc Package

AEU-7000/ Mont Blanc LED Package

AEU-7000/Impulsion Package

Handpieces



Aseptico Handpieces



20:1 Economy Dental Implant Handpiece (Max 70 Ncm)

Product Name	Part #
20:1 Reduction	AHP-85-SS



Mont Blanc 20:1 Reduction Contra Angle Handpiece

Product Name	Part #
20:1 Mont Blanc Reduction Contra Angle	AHP-85MB-X



20:1 Mont Blanc Handpiece w/Depth Gauge & Fiber Optics

Product Name	Part #
20:1 Mont Blanc w/Depth Gauge & Fiber Optics	AHP-85MBFO-CX



Straight AT1

Product Name	Part #
1:1 Straight Angled Handpiece	AT-1

Aseptico Accessory Products

General Dental Products



Autoclavable Replacement Tubing Set*

Product Name	Part #
Autoclavable Replacement Tubing Set	AE-23

*For the AEU-7000 and AEU-6000 series implant, surgery, & endo motors



Sodium Chloride - 0.9% 1000ml Bag

Product Name	Part #
Sodium Chloride	SODIUM CHLORIDE



Powered Trolley Cart

Product Name	Part #
Powered Trolley Cart	ATC-12v2

Terms & Conditions



Shipping, Return Policy, Implant Warranty, & Contact Information

How Can I Order Your Products?

Call our customer service department at 800-228-0477
Monday - Friday 8:00am to 5:00pm Mountain Time. You can
also go to www.ocobiomedical.com.

When Will My Shipment Arrive?

OCO Biomedical ships through FedEx exclusively, transit times
are as follows:

Overnight: Delivery next business day by 10:30am

2nd day: Delivery in 2 business days by 4:30pm

Ground: Delivery within 4-5 business days, by the end of day.

What is Your Return Policy?

All product must be in new and saleable condition (returned
in the original packaging) with the manufacturing seals intact.
Original invoice number and purchase date must accompany
product in order to process the return. Product returned
within 30 days of purchase can be returned for refund in the
form of original payment. Product returned 30-60 days after
purchase can be exchanged for same product or credit on
account only. Product returned 60-120 days after purchase
will require a 15% restocking fee and credit on account only.
Product returned 120 days-1 year after purchase will be
exchanged for same product (implants only). No products
can be returned after 1 year, with no exceptions. The returned
product must be enclosed in safe, protective packaging such
as a padded envelope or box. Please mark the outside of the
package to the attention of "Returns".

Return Address: OCO Biomedical Inc.,
Attn: Returns
9550 San Mateo Blvd. NE, Suite C, Albuquerque, NM 87113

Do You Have an Implant Warranty?

OCO Biomedical Dental Implants are warranted
unconditionally against defects in material and workmanship
for 5 years after production. The warranty covers only the
cost of the implant and does not cover professional fees or
restoration components. If a defect is found, please do not
use the item. Simply return the item after notifying OCO
Biomedical, and a replacement item will be delivered at no
additional charge. OCO Biomedical Dental Implants are
conditionally warranted if failure occurs within 2 years after
placement provided that: the implant was returned
un-cleansed, the implant was used with an OCO Biomedical
abutment, and the implant and abutment, if applicable, are
returned to OCO Biomedical accompanied by a completed
Implant Incident Report, including lot number, catalog
number and pre/post x-rays. The implant must be returned
within 30 days of removal/failure. A Implant Incident
Report may also be obtained through the Quality Control
Department.

Product Sterilization Disclaimer

Implants: OCO Biomedical opposes the resterilization
of implants after they have been removed from the
vial or placed in the jaw because of the possibility of
the transference of blood borne pathogens or other
contaminants.

All other surgical instrumentation: OCO Biomedical
recommends using a steam sterilization unit at 121
degrees Celsius for 15 minutes with a drying time of
30 minutes.

Contact Us

Phone

Toll Free: (800) 228-0477
International: +1 (505) 293-0025
Direct: (505) 293-0025
Fax: (505) 293-0447

Address

OCO Biomedical Inc.
9550 San Mateo Blvd. NE, Suite C
Albuquerque, NM, USA 87113

Hours

8:00 am - 5:00 pm MST
Monday - Friday

www.ocobiomedical.com